



### WEEKLY TIMER TEST

- Automatic weekly tests ensure that the EZ CARE + device is operating properly at all times.



### PROGRAMMING

- A simple online form is all that's needed to set up the station.
- No cumbersome on-device programming needed.



### MONITORED BATTERY LIFE

- When pressed, pendant will send battery life to base station. Will send an alert when the battery is low.
- 24 hour battery backup on base station. Will notify central when battery is low.



### OPTIONAL WIRELESS PERIPHERALS

- Smoke Alarm
- Carbon Monoxide Alarm
- Motion Detector
- These alarms and detectors provide accurate and vitally important inactivity alerts.



### CELLULAR CONNECTIVITY

- No phone line required
- Communicates over the GSM cellular network (T-Mobile or AT&T) using CID into a Sur-gard receiver.



### MECHANICS

- Dimensions: 8.12 in. x 5.5 in. x 1.25 in.
- Weight: Base unit 9.2 ounces, Pendant 0.4 ounces
- Operating Temperature: -50° C to 120° C
- Humidity: 0% to 90% non-condensing



### RANGE

- The connecting pendant has a range of up to 400 feet of open air from the base station.
- Easily connects up to 50 devices to one base station.