

LifeBeacon™



FREEROAM TECHNOLOGY

- FreeRoam Technology allows the customer to use the LifeBeacon wherever you are. *

*Requires adequate cellular telephone signal.



TWO-WAY VOICE

- Talk and listen to medical monitoring operators directly from your LifeBeacon.



INSTALINK

- Press the button and summon help in seconds, not minutes.*

*Average connection time is 28 seconds vs. traditional PERS at over a minute.



GPS

- LifeBeacon communicates with satellites to locate you in an emergency and then informs responders of your location.



INSTALERT

- LifeBeacon notifies loved ones immediately of your need for help via text or email along with a map link of your location.



AUDIBLE NOTIFICATIONS

- LifeBeacon will let you know when you've activated it, when to recharge it and more.



EXTENDED BATTERY LIFE

- LifeBeacon will last up to 10 days before needing to be recharged.*

*Battery life fluctuates depending on how often the emergency button is pressed. On average, non-emergency situations battery last around 10 days.



SIMPLICITY

- Set up is easy. No installation. Easy online dealer portal.



WATER RESISTANT

- The LifeBeacon can be used in the shower or worn in the rain.

ELECTRICAL

- 800 mAh Battery
- Current 0.9mA on standby and 1.3mAh avg. active current.
- Charges through magnetic charging contacts on the back of the device.
- Anticipated battery life up to 10 days in standby mode.

MECHANICAL

- Dimensions: 62mm x 41mm x 15.5mm
- Weight: 37g (without lanyard or bracelet)
- Gross weight of one complete box: 241g

COMMUNICATION

- The LifeBeacon communicates CID over Fibro Protocol via IP to a Sur-gard receiver.
- All mapping is done in the LiveFree cloud and is immediately updated to the central station.